Aven Tools

Safety Procedures

Safety Operating Guidelines For:

- Microscopes
- Digital Microscopes
- Magnifying Lamps
- Small Precision Hand Tools

Read and understand all safety procedures before operating

Tools for
Advancing
Innovation™

Website: aventools.com

Email: info@aventools.com

Phone: (734) 973-0099

Table of Contents

English	2-5
Español	6-9
Français	10-13
Italiano	14-17
Deutsch	18-21

Standard safety and operating guidelines for a microscope camera

A microscope camera is a device that allows you to capture images or videos of microscopic objects or specimens. It can be attached to a microscope or used as a standalone device. To ensure the safety of yourself and others, and to prevent damage to the equipment, please follow these standard safety and operating guidelines when using a microscope camera:

- Before using the microscope camera, read the user manual carefully and familiarize yourself with its features and functions.
- Do not use the microscope camera in wet, dusty, or humid environments, or near sources of heat or electromagnetic interference.
- Do not expose the microscope camera to direct sunlight or strong light sources, as this may damage the sensor or affect the image quality.
- Do not touch the lens or sensor of the microscope camera with your fingers or any sharp or abrasive objects, as this may scratch or damage them.
- Do not apply excessive force or pressure to the microscope camera or its cables, as this may cause malfunction or breakage.
- Do not disassemble or modify the microscope camera or its accessories, as this may void the warranty and cause injury or damage.
- If the microscope camera is attached to a microscope, make sure it is securely mounted and aligned with the eyepiece or objective lens of the microscope.
- If the microscope camera is used as a standalone device, make sure it is placed on a stable and level surface, and adjust the focus and zoom knobs until you get a clear image of the object or specimen.
- To capture images or videos with the microscope camera, connect it to a computer or monitor using the appropriate cable or wireless connection, and use the software provided by the manufacturer or a compatible software application.
- To store or transfer images or videos from the microscope camera, use a memory card, USB flash drive, external hard drive, cloud service, or other suitable storage device or method.
- To clean the microscope camera, use a soft cloth or lens tissue moistened with water or alcohol. Do not use any solvents or chemicals that may damage the surface or components of the microscope camera.
- To store the microscope camera when not in use, place it in a protective case or bag, and keep it in a cool, dry, and dark place. Avoid extreme temperatures and humidity levels that may affect its performance or lifespan.

By following these standard safety and operating guidelines for a microscope camera, you can ensure your safety and satisfaction while using this device. If you have any questions or problems with the microscope camera, please contact the manufacturer or an authorized service center for assistance.

Standard safety operating guidelines for magnifying lamps

Magnifying lamps are devices that use a lens and a light source to magnify and illuminate small objects or details. They are commonly used for various purposes such as reading, crafting, inspecting, soldering, and medical examinations. However, magnifying lamps also pose some potential hazards if not used properly. Therefore, it is important to follow these standard safety operating guidelines when using magnifying lamps:

- Do not modify or alter the device in any way.
- Always plug the magnifying lamp into a grounded outlet with the correct voltage and current rating. Do not use extension cords or adapters that are not compatible with the device.
- Do not place the magnifying lamp near flammable materials, liquids, or gases. Keep it away from heat sources, direct sunlight, and moisture. Do not expose it to extreme temperatures or humidity.
- Do not touch the lens, the light bulb, or any metal parts of the magnifying lamp when it is on or shortly after turning it off. They may be hot and cause burns or electric shocks.
- Do not look directly at the light source of the magnifying lamp or use it to view bright objects such as the sun or lasers. This may damage your eyes or cause temporary blindness.
- Adjust the position and angle of the magnifying lamp carefully and securely. Do not force or bend the arm or the joints of the device. Do not move or carry the magnifying lamp by holding the lens or the light bulb.
- Clean the lens regularly with a soft cloth. Do not use abrasive materials, solvents, or alcohol to clean them. Do not immerse the device in water or any other liquid.
- Replace the light bulb when it is dim, flickering, or burned out. Use only the same type and wattage of bulb as specified by the manufacturer. Turn off and unplug the device before replacing the bulb. Wear gloves and handle the bulb carefully to avoid breaking it.
- If you notice any damage, malfunction, or abnormality with the magnifying lamp, stop using it immediately and contact the manufacturer or a qualified service technician for repair or replacement. Do not attempt to fix it yourself unless you are authorized and trained to do so.

By following these standard safety operating guidelines for magnifying lamps, you can ensure a safe and effective use of this device for your intended purposes.

Standard Safety Operating Guidelines for Microscopes/Video Inspections Systems

Microscopes are valuable tools for scientific research, manufacturing and education, but they also require proper handling and maintenance to ensure their optimal performance and safety. The following guidelines provide some basic tips on how to use and care for microscopes in a professional manner.

- Before using a microscope, make sure it is plugged into a suitable power source and that the power cord is not damaged or frayed. Do not attempt to disassemble or modify the microscope, as this may cause electric shock or damage to the instrument.
- Adjust the microscope to suit your posture and comfort level. Use ergonomic accessories such as extended eye tubes, optical wedges, flexible adapters, or adjustable stands if available. Avoid bending your neck, back, or wrists excessively while observing specimens. Take frequent breaks to relax your eyes and muscles.
- Clean the lenses with lens paper or a soft cloth moistened with a small amount of alcohol-ether mixture or isopropyl alcohol. Do not use water, solvents, or abrasive materials, as they may scratch or damage the lenses. Do not touch the lenses with bare fingers, as this may leave fingerprints or oils that affect the image quality.
- Clean the microscope body and stand with a moist, soft cloth with a small amount of detergent. Do not use water on plastic surfaces, as this may cause cracking or discoloration. Do not use paint thinner or other solvents, as they may damage the coating or parts of the microscope. Dry all surfaces after cleaning.
- Always cover the microscope with a dust cover when not in use. Store it in a dry place away from direct sunlight, heat sources, or moisture. In humid or moist environments, it is advisable to store the microscope in a waterproof container with a drying agent.

Standard safety and operating guidelines for small precision hand tools

Small precision hand tools are instruments that are used to perform delicate tasks that require high accuracy and skill. These tools include screwdrivers, pliers, tweezers, cutters, scissors, knives, files, drills, hammers, and saws. They are often used in fields such as electronics, jewelry, watchmaking, woodworking, and crafts.

To use these tools safely and effectively, it is important to follow some standard guidelines:

- Before using any tool, inspect it for any damage or defects. Do not use a tool that is broken, dull, rusty, or dirty.
- Choose the right tool for the job. Do not use a tool that is too big or too small for the task. Do not use a tool for a purpose that it is not designed for.
- Wear appropriate personal protective equipment (PPE) such as gloves, goggles, masks, and aprons. PPE can protect you from injuries such as cuts, burns, electric shocks, and eye damage.

- Work in a well-lit and ventilated area. Avoid working near flammable or explosive materials. Keep your work area clean and organized. Store your tools in a safe and dry place when not in use.

- Handle your tools with care and respect. Do not force or over-tighten a tool. Do not drop or throw a tool. Do not leave a tool unattended or in a position where it can fall or cause harm.

- Follow the manufacturer's instructions and recommendations for using and maintaining your tools. Sharpen your tools regularly and lubricate them if needed. Clean your tools after each use and disinfect them if they come in contact with blood or other bodily fluids.

- Seek medical attention immediately if you suffer an injury from using a tool. Report any accidents or incidents to your supervisor or manager.

By following these standard safety and operating guidelines for small precision hand tools, you can ensure that you perform your work efficiently and safely.

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